Brief Introduction of Y41SE STR Detection Kit



1. Loci Arrangement Diagram of Y41SE STR Detection Kit





2. Features of Y41SE STR Detection Kit

Six-color fluorescence detection platform, can amplify 40 Y-STR loci including 33 low-mutant loci (mutant rate less than 1%) and 7 rapid mutant loci (mutant rate higher than 1%) simultaneously, and one Y indel with a very low mutant rate.

2. Y-STR database construction, suitable for male family investigation, featuring a strong individual recognition ability, and having all 35 loci required and preferred in the Y library standard.

3. Y41SE contains **12 mini-STR** (amplifying length under 220bp), which can acquire more information from DNA degraded samples.

4. **Strong compatibility,** compatible with the direct amplification of various samples and the amplification of purified DNA.

5. Rapid amplification, the whole process is 70min theoretically.



4.1 Components and Contents of Y41SE STR Detection Kit

Components of PCR Amplification Kit	100 Reactions (Y series)
2.5× PCR Mix	1000μL× 1
5× Primers	500μL× 1
PCR Grade Water	1000µL× 1
Control DNA 9948(0.5ng/µL)	20μL× 1
Taq DNA Polymerase	100μL× 1

Components of PCR Detection Kit	100 Reactions (Y series)	
Allelic Ladder	25μL× 1	
Orange-600	120μL× 1	
Components of Matrix Kit	Free of charge for the first time	
6 Dye Matrix	50μL× 1	



4.2 Reaction System of Y41SE STR Detection Kit

Component	25µL System	10µL System
2.5× PCR Mix	10.0	4.0
5×Primers	5.0	2.0
PCR Grade Water	8.0	2.6
Control DNA 9948(0.5ng/µL)	1.0	1.0
Taq DNA Polymerase	1.0	0.4
Total Volume	25.0	10.0

Note: Customers are not recommended to adopt reaction systems under 10µl.



4.3 Reaction System and Procedures of Y41SE STR Detection Kit







WE ARE JUST ON THE WAY THANK YOU.

Contact Us Website: http://www.dnatechlab.com/



江苏苏博生物医学科技南京有限公司